


SOP for Formlabs Custom trays and Wax Rims


By:
Katherine Lala

Nesting: software is PreForm (blue butterfly)

- when the dialog box opens pick the printer of choice
- resin is Custom Tray
- layer thickness set to fastest print 200 microns
- apply
- drag and drop files from the windows directory \dima Print Cast
- Orient print with the intaglio surface away from the plate (so that minimal supports get on it.)



- if there are a fewer than five prints lay them flat to the plate as above



- if there are more than that, stand them nearly vertical to the plate distal side to plate (with a slight tilt with the intaglio side away from the plate) as above
- Supports:
 - raft type = Mini rafts
 - no raft labels
 - density = 1.00
 - touchpoint size = 0.30mm
 - internal supports off
 - Auto-Generate All
- **Check for and supports on the intaglio surfaces.** Remove if necessary.
 - For easier visibility select the model. It will turn blue and the supports will be green.
- Layout All
- Check the slices.
- Save file as month-day-70(2/3)-#-initials
- Check printability
- use the Orange Printer button to send to the printer.

Printing:

- Remove the resin tank and while holding upside down give it a left right shake. Put it back.
- Open the vent on resin tank
- check the plate (clean and proper position)
- in the tiger mix the resin
- verify that you have the correct file and hit print on the printer. Follow prompts.

Removal and Curing:

- remove plate
- at the processing station use the super bendy black handle scraper to remove the prints
- place them in the dirty IPA for 280 seconds in the ultrasonic
- blow out between supports and print

- place in clean IPA for 280 seconds
- give a light blow off and air dry for 30 minutes
- 30 minutes in the Accu-Cure
- remove supports and clean up with clippers

formlabs reference <https://dental.formlabs.com/indications/custom-trays-guide/>